

21 October 2026

Tenofovir Disoproxil Viatriis (tenofovir disoproxil maleate) 300 mg film coated tablet

Dear Healthcare Professional,

We are writing to inform you of an important update regarding Tenofovir Disoproxil Viatriis (tenofovir disoproxil maleate) 300 mg film-coated tablets, sponsored by Viatriis Ltd.

To support continuity of supply, Medsafe has granted consent for Tenofovir Disoproxil Viatriis (tenofovir disoproxil maleate) 300 mg film-coated tablets manufactured for the UK market to be temporarily supplied in New Zealand.

The UK and New Zealand products are identical in formulation. The only difference between the UK and New Zealand presentations is a slight difference in the tablet debossing (see table below). Viatriis intends to distribute approximately two months' supply of the UK products labelled with NZ compliant labels to support continuity of supply in New Zealand, during which time both the UK and New Zealand presentations will be supplied concurrently.

The product label for the UK product will be identical to the registered New Zealand product label.

The batch number of the UK product is: 8222253.

Adverse Event Reporting & Medical Enquiries

Please report any suspected adverse events via email to Medinfo_anz@viatriis.com. Alternatively, suspected adverse events may be reported to the Centre for Adverse Reactions Monitoring (CARM) online at <https://pophealth.my.site.com/carmreportnz/s/> reporting or by email to CARMreport@health.govt.nz. Please direct any medical enquiries to Viatriis or report any suspected adverse drug reactions to Viatriis via telephone on 0800 168 169 or by email at medinfo_anz@viatriis.com.

Yours sincerely,

Elissa Soo

Elissa Soo

Senior Medical Affairs Specialist

Overview of the difference between the New Zealand and UK Tenofovir Disoproxil Viatris

New Zealand tablet appearance	The UK tablet appearance
Light blue coloured, round, biconvex, film coated tablets debossed with 'TM300' on one side of the tablet and 'M' on other side. Dimensions (diameter): 12.20 mm $\pm$ 0.20 mm	Light blue coloured, round, biconvex, film coated tablets debossed with 'TN245' on one side of the tablet and 'M' on other side. Dimensions (diameter): 12.20 mm $\pm$ 0.20 mm